## INSTRUCTIONS FOR USE

**CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.**

### Device Description

The ATLAS<sup>®</sup> GOLD PTA Dilatation Catheter is a high performance balloon catheter consisting of an over-the-wire catheter with a balloon fixed at the distal tip. The proprietary non-compliant, low profile balloon is designed to provide consistent balloon diameters and lengths even at high pressures. Two radiopaque markers delineate the working length of the balloon and aid in balloon placement. The coaxial catheter includes a tapered atraumatic tip to facilitate advancement of the catheter to and through the stenosis. The proximal portion of the catheter includes a female luer lock hub connected to the inflation lumen, and a female luer-lock hub connected to the guidewire lumen. The over-the-wire catheter is compatible with .035" guidewire and is available in 80cm and 120cm working lengths.

Packaged with every product is a profile reducing sheath that is positioned over the balloon for protection before use. A re-wrapping tool is also provided on the catheter shaft to aid in re-wrap/refolding of the balloon.

This product is not manufactured with any natural rubber latex.

### Indications for Use

ATLAS<sup>®</sup> GOLD PTA Dilatation Catheter is indicated for use in Percutaneous Transluminal Angioplasty of the peripheral vasculature, including the iliac arteries and iliac and femoral veins, and for the treatment of obstructive lesions of native or synthetic arteriovenous dialysis fistulae. This device is also indicated for post-dilatation of stents and stent grafts in the peripheral vasculature. This catheter is not for use in coronary arteries.

### Contraindications

None known.

### Warnings

1. Contents supplied **STERILE** using ethylene oxide (EO). Non-Pyrogenic. Do not use if sterile barrier is opened or damaged. Single patient use only. Do not reuse, reprocess, or re-sterilize.
2. This device has been designed for single use only. Reusing this medical device bears the risk of cross-patient contamination as medical devices – particularly those with long and small lumina, joints, and/or crevices between components – are difficult or impossible to clean once body fluids or tissues with potential pyrogenic or microbial contamination have had contact with the medical device for an indeterminable period of time. The residue of biological material can promote the contamination of the device with pyrogens or microorganisms which may lead to infectious complications.
3. Do not resterilize. After reesterilization, the sterility of the product is not guaranteed because of an indeterminable degree of potential pyrogenic or microbial contamination which may lead to infectious complications. Cleaning, reprocessing, and/or reesterilization of the present medical device increases the probability that the device will malfunction due to potential adverse effects on components that are influenced by thermal and/or mechanical changes.
4. To reduce the potential for vessel damage, the inflated diameter and length of the balloon should approximate the diameter and length of the vessel just proximal and distal to the stenosis.
5. To reduce the potential for stent or stent graft damage and/or vessel damage from the stent or stent graft, the diameter of the balloon should be no greater than the diameter of the stent or stent graft. Refer to the stent or stent graft IFU for safety information including the WARNINGS, PRECAUTIONS, and potential ADVERSE EFFECTS regarding the use of balloon post-dilatation.
6. When the catheter is exposed to the vascular system, it should be manipulated while under high-quality fluoroscopic observation. Do not advance or retract the catheter unless the balloon is fully deflated. If resistance is met during manipulation, determine the cause of the resistance before proceeding. Applying excessive force to the catheter can result in tip breakage or balloon separation.
7. Do not exceed the RBP recommended for this device. Balloon rupture may occur if the RBP rating is exceeded. To prevent over pressurization,

use of a pressure monitoring device is recommended.

8. After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state, and federal laws and regulations.

### Precautions

1. Carefully inspect the catheter prior to use to verify that catheter has not been damaged during shipment and that its size, shape, and condition are suitable for the procedure for which it is to be used. Do not use if product damage is evident.
2. The ATLAS<sup>®</sup> GOLD Catheter shall only be used by physicians trained in the performance of Percutaneous Transluminal Angioplasty.
3. The minimal acceptable sheath French size is printed on the package label. Do not attempt to pass the PTA catheter through a smaller size sheath introducer than indicated on the label.
4. Do not remove the guidewire in situ to shoot contrast through the wire lumen or perform a wire exchange. If the wire is removed while the balloon catheter is situated in tortuous anatomy, the risk of kinking the catheter is increased.
5. Use the recommended balloon inflation medium (a range of 30-50% contrast medium / a range of 50-70% sterile saline solution). It has been shown that a 30/70% contrast/saline ratio has yielded faster balloon inflation / deflation times.
6. Never use air or other gaseous medium to inflate the balloon.
7. If resistance is felt during post procedure withdrawal of the catheter through the introducer sheath, determine if contrast is trapped in the balloon with fluoroscopy. If contrast is present, push the balloon out of the sheath and then completely evacuate the contrast before proceeding to withdraw the balloon.
8. If resistance is still felt during post procedure withdrawal of the catheter, it is recommended to remove the balloon catheter and guidewire/introducer sheath as a single unit.
9. Do not continue to use the balloon catheter if the shaft has been bent or kinked.
10. Prior to re-insertion through the introducer sheath, the balloon should be wiped clean with gauze, rinsed with sterile normal saline, and refolded with the balloon re-wrap tool. Balloon re-wrapping should only occur while the balloon catheter is supported with a guidewire.

### Potential Adverse Reactions

The complications which may result from a peripheral balloon dilatation procedure include:

- Acute thrombotic occlusion
- Additional intervention
- Allergic reaction to drugs or contrast medium
- Aneurysm or pseudoaneurysm
- Arrhythmias
- Balloon rupture
- Balloon getting stuck on stent
- Distal embolization (PE)
- Hematoma
- Hemorrhage, including bleeding at the puncture site
- Hypotension/hypertension
- Inflammation
- Leg edema
- Occlusion
- Pain or tenderness
- Pneumothorax or hemothorax
- Sepsis/infection
- Shock
- Short term hemodynamic deterioration
- Stent disruption or dislodgement with balloon insertion
- Stroke
- Thrombosis
- Vessel dissection, perforation, rupture, or spasm

### Directions for Use

#### Handling & Storage

Store in a cool, dry, dark place. Do not store near radiation or ultra-violet light sources.

Rotate inventory so that the catheters and other dated products are used prior to the "Use By" date.

Do not use if packaging is damaged or opened.

#### Equipment for Use

- Contrast medium
- Sterile saline solution
- Luer lock syringe/inflation device with manometer (10 ml or larger)
- Appropriate introducer sheath and dilator set
- .035" guidewire

#### Dilatation Catheter Preparation

1. Remove Catheter from package.

2. Verify the balloon size is suitable for the procedure and the selected accessories accommodate the catheter as labeled.
3. Remove the balloon guard by grasping the balloon catheter just proximal to the balloon and with the other hand, gently grasp the balloon protector and slide distally off of the balloon catheter.
4. Slide the re-wrap tool to the proximal end of the catheter shaft.
5. Prior to use, the air in the balloon catheter should be removed. To facilitate purging, select a syringe or inflation device with a 10 ml or larger capacity and fill approximately half of it with the appropriate balloon inflation medium (30-50% contrast medium / 50-70% sterile saline solution). Do not use air or any gaseous medium to inflate the balloon.
6. Connect a stopcock to the balloon inflation female luer hub on the dilatation catheter.
7. Connect the syringe to the stopcock.
8. Hold the syringe with the nozzle pointing downward, open the stopcock, and aspirate for approximately 15 seconds. Release the plunger.
9. Repeat step #8 two more times or until bubbles no longer appear during aspiration (negative pressure). Once completed, evacuate all air from the barrel of the syringe/inflation device.
10. Prepare the wire lumen of the catheter by attaching a syringe to the wire lumen hub and flushing with sterile saline solution.

### Use of the ATLAS<sup>®</sup> GOLD PTA Dilatation Catheter

1. Backload the distal tip of the ATLAS<sup>®</sup> GOLD Catheter over the pre-positioned guidewire and advance the tip to the introduction site.
2. Advance the catheter through the introducer sheath and over-the-wire to the site of inflation. If the stenosis cannot be crossed with the desired dilatation catheter, use a smaller diameter catheter to pre-dilate the lesion to facilitate passage of a more appropriately sized dilatation catheter.
3. Position the balloon relative to the lesion to be dilated, ensure the guidewire is in place and inflate the balloon to the appropriate pressure.
4. Apply negative pressure to fully evacuate fluid from the balloon. Confirm that the balloon is fully deflated under fluoroscopy.
5. While maintaining negative pressure and the position of the guidewire, withdraw the deflated dilatation catheter over-the-wire through the introducer sheath. Use of a gentle clockwise motion may be used to help facilitate catheter removal through the introducer sheath.

### Balloon Reinsertion

**Precaution:** Do not continue to use the balloon catheter if the shaft has been bent or kinked.

**Precaution:** Prior to re-insertion through the introducer sheath, the balloon should be wiped clean with gauze, rinsed with sterile normal saline, and refolded with the balloon re-wrap tool. Balloon re-wrapping should only occur while the balloon catheter is supported with a guidewire.

1. Load the balloon catheter onto a guidewire.
2. Advance the balloon re-wrap tool over the catheter to the proximal end of the balloon.
3. Grasp the catheter shaft just proximal to the balloon with one hand, and with the other hand gently slide the rewrap tool over the balloon to the catheter tip and then back over the balloon to the catheter.
4. Slide the re-wrap tool to the proximal end of the catheter shaft.
5. Advance the balloon catheter over the pre-positioned guidewire to the introduction site and through the introducer sheath. If resistance is encountered, replace the previously used balloon catheter with a new balloon.
6. Continue the procedure according to the "Use of the ATLAS<sup>®</sup> GOLD PTA Dilatation Catheter" section herein.

**Warning:** After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practices and applicable local, state and federal laws and regulations.

### Clinical Experience

#### Design

A physician sponsored, single center, retrospective study was conducted to assess the safety and effectiveness of the ATLAS<sup>®</sup> GOLD PTA Dilatation Catheter for balloon dilatation in iliac and femoral deep venous procedures up to one year post-procedure. Sixty-one patients treated with the ATLAS<sup>®</sup> GOLD PTA Dilatation Catheter of the iliofemoral veins were included.

The primary safety endpoint was defined as intra-procedural freedom from major adverse events (MAE) of  $\geq 95\%$ , defined as acute thrombosis, perforation, or device-related complications. The secondary endpoints were in-hospital and 1 month freedom from MAE; minimal luminal area (MLA) post-stent minus MLA pretreatment at compression site using intravascular ultrasound (IVUS)

[absolute MLA gain]; clinical improvement at one year (when available) +/- 2 months, defined as reported by patients (less swelling and/or pain in the treated limb); and patency by duplex ultrasound at one month and 1 year when data was available.

**Patient Demographics**

The mean age of patients was 61.9 years, and 60.7% were females. Eighteen percent of patients had a history of deep vein thrombosis (DVT) and 13.1% of patients had history of pulmonary embolism (PE). The majority of patients did not have a history of smoking, diabetes, or peripheral arterial disease, nor presented with a hypercoagulable state or deep vein reflux. A majority of patients presented with heaviness or pain and swelling of the leg. One patient had an IVC filter.

IVUS Common Iliac Vein Area Ipsilateral Distal Reference (mm <sup>2</sup> )	40	226.7	85.0-455.0
IVUS Post-Stent Area at Confluence (mm <sup>2</sup> )	39	218.8	87.0-433.0

**Conclusion**

The ATLAS<sup>®</sup> GOLD PTA Dilatation Catheter met 100% freedom from intraprocedural adverse events related to the study device, and exceeded the 95% benchmark for the primary safety endpoint. No balloon ruptures, vessel perforations, or balloon-related intravascular events occurred. At one month and one year there were no adverse events related to ATLAS<sup>®</sup> GOLD PTA Dilatation Catheter use, with marked symptom improvement in the majority of patients. Additionally, imaging analysis showed that stent dilatation yielded excellent stent expansion at the level of compression and the MLA obtained was nearly equivalent to the reference diameter.

**Warranty**

Bard Peripheral Vascular warrants to the first purchaser of this product, that this product will be free from defects in materials and workmanship for a period of one year from the date of first purchase and liability under this limited product warranty will be limited, to repair or replacement of the defective product, in Bard Peripheral Vascular's sole discretion, or refunding your net price paid. Wear and tear from normal use or defects resulting from misuse of this product are not covered by this limited warranty.

**TO THE EXTENT ALLOWABLE BY APPLICABLE LAW, THIS LIMITED PRODUCT WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT WILL BARD PERIPHERAL VASCULAR BE LIABLE TO YOU FOR ANY INDIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES RESULTING FROM YOUR HANDLING OR USE OF THIS PRODUCT.**

Some countries do not allow an exclusion of implied warranties, incidental or consequential damages. You may be entitled to additional remedies under the laws of your country.

**An issue or revision date and revision number for these instructions are included for the user's information on the last page of this booklet.**

In the event 36 months have elapsed between this date and product use, the user should contact Bard Peripheral Vascular to see if additional product information is available.

Patient Demographics				
		N (N=61)	Mean	Range
Age	Female	37	60.9	23.0-85.0
	Male	24	63.3	24.0-87.0
	ALL	61	61.9	23.0-87.0
BMI (kg/m2)	Female	37	32.0	19.1-45.3
	Male	24	33.5	20.5-49.2
	ALL	61	32.6	19.1-49.2

Variable		N (N=61)	%
History of DVT	No	50	82.0%
	Yes	11	18.0%
Variable		N (N=11)	%
Location of DVT	*	4	36.4%
	Left popliteal	2	18.2%
	Bilateral Unnamed	1	9.1%
	Left Axillary	1	9.1%
	Left Lower Extremity	3	27.3%
Variable		N (N=61)	%
History of PE	No	53	86.9%
	Yes	8	13.1%
Hypertension	No	26	42.6%
	Yes	35	57.4%
Smoking	No	40	65.6%
	Yes	21	34.4%
Diabetes Mellitus	No	52	85.2%
	Yes	9	14.8%
Peripheral Arterial Disease	No	58	95.1%
	Yes	3	4.9%
Inferior Vena Cava Filter	No	60	98.4%
	Yes	1	1.6%
Deep Vein Reflux	No	43	70.5%
	Yes	9	14.8%
	*	9	14.8%
Heaviness/Pain in Leg	No	18	29.5%
	Unilateral	43	70.5%
	Bilateral	0	0.0%
Swelling in Lower Leg	No	14	23.0%
	Unilateral	47	77.0%
	Bilateral	0	0.0%
Hypercoagulation State	No	60	98.4%
	Yes	1	1.6%

\*not recorded

**Procedural Variables**

Balloons used for pre-dilatation ranged from 12-14 mm in diameter and were inflated to 3-14 atm. Balloons used for post-dilatation ranged from 12-18 mm in diameter and were inflated to 2-20 atm. All patients received a stent.

**Results**

There were no intraprocedural complications for 96.8% of patients treated with ATLAS<sup>®</sup> GOLD PTA Dilatation Catheter. One patient experienced stent migration during stent deployment, which was unrelated to the use of ATLAS<sup>®</sup> GOLD PTA Dilatation Catheter. The IVUS MLA at the confluence post-stenting was markedly improved.

Intraprocedural Complications		
Variable	N (N=61)	%
No intraprocedural complications	60	96.8%
Distal embolization (causing documented PE)	0	0%
Acute thrombotic occlusion	0	0%
Balloon rupture	0	0%
Balloon getting stuck on stent	0	0%
Stent disruption or dislodgement with balloon insertion	0	0%
Stent disruption not related to study device	1	1.6%

Intraprocedural IVUS Measurements			
Variable	N (patients with IVUS recorded)	Mean	Range
MLA at compression (mm)	40	72.0	27.0-174.0
Percent compression by IVUS	40	66.2	15.9-91.1

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**Band and Atlas are trademarks and/or registered trademarks of C.R. Bard, Inc. or an affiliate.**
Bard and Atlas are C.R. Bard, Inc. or its affiliated company's trademarks and/or registered trademarks.
**Copyright © 2019 C. R. Bard, Inc. All rights reserved. Printed in the USA.**
Copyright © 2019 C. R. Bard, Inc. Tous droits réservés. Imprimé aux États-Unis.
Copyright © 2019 C. R. Bard, Inc. Alle Rechte vorbehalten. Gedruckt in den USA.
Copyright © 2019 C. R. Bard, Inc. Tutti i diritti riservati. Stampato negli Stati Uniti.